



## MATa

Saturated culture in 200µl YPD  
+ cloNAT in 96-well plate




## MATα


Saturated culture in 35ml YPD  
+ hygromycin in 50 ml tube



**MATING:** Omnitray with  
35ml of solid YPD





## MATa/α



**DIPLOID SELECTION:**  
Omnitray with 35ml of solid  
YPD + cloNAT + hygromycin

1536-position  
pin tool



**SCREEN OF PPIs:**  
Omnitray with 35ml of solid  
synthetic media -ade +MTX