



S3 Fig. Mouse RNase L but not human RNase L binds to L*.

Surface plasmon resonance (SPR) experiments were performed using Biacore technology with purified recombinant L*, mouse RNase L and human RNase L.

A. His6-L* was immobilized on Ni-NTA chips and human RNase L was injected at the indicated concentrations.

B. His6-L* was immobilized on Ni-NTA chips and mouse RNase L was injected at the indicated concentrations.

C. GST-mouse RNase L was immobilized on anti-GST CM5 sensor chip and L* was injected at the indicated concentrations.

RU: response units